

CHEMOTHERAPY / THERAPEUTIC CHAIRS

Genesys Chemotherapy/Therapeutic chairs are designed focusing on Quality, Comfort, Safety and Hygiene. The stylish design comes with special features that distinguish it as a leading oncology chair. The remote control handset allows the chair to be positioned accordingly with little effort required by both operator and patient. Quick conversion to shock position within 10 seconds of actuation. Designed with two actuators (three motors are optional) for the easy movement of backrest & leg rest. Automatic synchronization between the back & leg supports, which increases the donor comfort while flattening. The Arm rest can be optimally set at the desired height and also can be moved to any direction. Designed to shift easily to Trendelberg position if any reactions are observed on the patient, then can be converted back to flat bed position. The height adjustment is done vertically with the help of reliable heavy duty lifting columns, thereby not causing any change in the profile of the donor. Logically placed buttons on remote enables easy to operate quick adjustments based on individual preferences. Integrated rotatable IV pole & Trays for placing various items. Large soft type wheels allow easy movement in the desired direction and facilitate easy cleaning without having to shift the equipment.

TECHNICAL SPECIFICATIONS

- ▣ Fold back upholstered armrests
- ▣ Lifting capacity : Up to 200Kgs
- ▣ Non-abrasive and non-deformable soft PU foam
- ▣ Extremely adjustable including horizontal function
- ▣ Excellent clinical access to assist cannulation and other procedures
- ▣ Adjustable headrest, Cushion body for ultimate comfort
- ▣ Adjustable foot rest with removable transparent plastic cover
- ▣ Removable and height adjustable pillow
- ▣ Stylish design facilitates easy cleaning
- ▣ Automatic synchronization between the back & leg supports
- ▣ Initiation of motor movement by wired remote with long lifetime- Up to 4 easy-touch buttons- Protection class: IP66 - 1m straight cable + 20 cm coiled cable- Hook for easy access
- ▣ Telescopic IV stand adjustable both sides with trays to hold BP Apparatus, IV bottles etc.
- ▣ Tissue / Towel holder is provided
- ▣ Compact, robust and reliable design
- ▣ Lifting capacity up to 200 kgs
- ▣ Construction: Mild steel with Epoxy powder coating
- ▣ Lockable soft type castor wheels
- ▣ Power : 220V/50Hz

