

# GROSSING WORKSTATION

"Genesys" Grossing workstation is totally of SS 304 Grade Stainless Steel construction that provides a more durable and long-lasting finish. All workstations are factory-fitted with exhaust system that can be connected to any in-house ventilation system. Our unique Back-Draft ventilation system removes hazardous fumes, vapors and odors from the work area, so you can be assured you're working in a safe environment.

The workstation is made of stainless steel from top-to-bottom that will withstand the harshest laboratory disinfectants. These models are self standing and sturdy unit. This unit maintains the air balancing by Negative Pressure with a suction velocity of 100 FPM with low noise level of 65 dB "A" Scale.

## TECHNICAL SPECIFICATIONS

The special features will be provided for the Grossing Workstation and also can be further customized to meet your needs by adding any of the following options;

- o Illumination - 4 Nos. of LED Lights,
- o 2 Nos. of universal 6/16A additional power sockets
- o Pull-out Writing Ledge below the table
- o Pull-out Pencil Drawer below the table
- o Paper Towel Holder inside the workstation
- o Shelf and Adjustable Bracket inside the work area
- o Exhaust Air Grill on the table and just below the storage platform
- o Storage Plat form is 8" width at inside the work area
- o Hand Water Spray Rinse on the table
- o Magnetic Tool Bar at rear panel
- o Grossing Perforated Tray on the Sink
- o Small size Polypropylene Platform of white colour for grossing process
- o Floor leveling screws at bottom of the stand
- o Supporting stand for weighing balance. *(optional)*
- o Monitor bracket designed to hold a flat-screen display monitor. *(optional)*
- o Keyboard holding bracket *(optional)*
- o CPU Mounting bracket. *(optional)*
- o 1/2 hp commercial waste disposal system.
- o Foot operated faucet control for water. *(optional)*
- o Halogen light mounted on a flexible arm. *(optional)*
- o Camera mount facilities like holder and splash shield. *(optional)*
- o Adjustable Magnifying lamp fixture. *(optional)*
- o Hands free soap dispenser. *(optional)*
- o Digital weighing scale and Centimeter ruler. *(optional)*
- o Computer system : All in one Computer system for Grossing Workstation, Wall Mountable, Intel core processor, Graphics / DVDRW / Wireless Keyboard & Mouse. *(optional)*
- o HD Camera, with feasibility of Ethernet, USB & Comport to capture the video / image by using the camera when specimen is being on processing in grossing workstation. *(optional)*

## MODEL SELECTION

Catalog No.	Working space Dimensions (WDH) in Inches	Suction Velocity (FPM)
GWS - 500H	60 x 24 x 30	110
GWS - 700	45 x 45 x 85	175



### MATERIAL OF CONSTRUCTION FOR THE GROSSING WORKSTATION:

**Main Hood & Bottom Support:** The complete cabinet is made out of thick heavy gauge SS 304 Grade Stainless Steel with matt finish. The necessary SS 304 heavy duty tubular supports and sides screen will be provided at the top. Cupboard is provided at the bottom side of the unit.

**Working Surface:** The work surface will be covered with thick heavy gauge SS 304 Grade Stainless Steel with Matt finish. Work surfaces are recessed and sloped to facilitate drainage into the extra large integral sink and also the Grossing Perforated Tray will be provided on the sink

## OPTIONAL ACCESSORIES

- o Perimeter rinse for free flow of water during grossing section
- o Formalin Monitor for assessing the level of formalin use in the room
- o Safety splash shield
- o Leg Frame
- o Eyewash assembly
- o Formalin Carboy
- o Glove box holder
- o IQ, OQ, PQ documentation